



suremark[®]

Precision Skin Marking Systems 

SINCE 1994
nltec

Celebrating 25 years as a part of the care-giving community.



We at Suremark® are grateful for your service to those in need of healing and positive energy. With this in mind, we have focused intensely on creating new **PRECISION SKIN MARKING PRODUCTS** that not only address the needs of technologists and doctors, but also provide a few moments of levity from the stress of medical procedures. With the overwhelming success of our ComfortMark™ and BeautyLine™ products, we have seen what many health care professionals are calling the "Suremark® Butterfly Effect." We've carried this good feeling into our newest product offerings for MRI, radiation therapy, and children's CT scan procedures.

The butterfly is one of the most emblematic symbols of personal transformation. In many traditions throughout the world, it represents renewal, rebirth, health and healing. Secondary meanings include lightness of being, playfulness, joy, and tuning into the emotional or spiritual. For Suremark®, the butterfly represents caring and understanding; that patients and their care-givers find themselves in difficult situations, and a kind word with a beautiful butterfly can go a long way.

Suremark® was the first to develop and offer a complete line of environmentally safe, **NON-METALLIC** precision skin markers to the medical industry, and our **LATEX-FREE** labels are designed and manufactured to be the finest quality available anywhere. From our advanced label adhesive, to our precise marker placement, Suremark® labels are guaranteed to be the best performing markers you have ever used. You can be confident that all Suremark® marking systems will stay where you put them, the first time, and show up clearly and consistently in all of your medical imaging applications.

Through our experience, Suremark® has implemented many manufacturing techniques that greatly improve efficiency, resulting in substantial savings over competing brands, and all Suremark® labels are backed by an unconditional money-back guarantee within one year of purchase.



Featured Suremark® Products



MRI Liquimark™: Our Liquimark MRI markers are NON-METALLIC and precision designed to the specific needs of MRI technologists and associated medical professionals.



Vision Line™ Kids: Give kids a choice when they feel like they don't have one with these fun, lighthearted precision skin markers. Puppies, kittens, pandas, and more for kids to choose their favorite!



Theramark™: Designed specifically for use on radiation therapy immobilization masks, the Suremark® Theramark™ is a precision label that allows technologists to accurately mark specific points for radiation treatment.



Scaneez® Cone Beam Markers: Stop wasting time and money duplicating dentures with radiopaque acrylics, metal bb's, gutta percha, etc. Scaneez® markers are designed to save time and money!



Tomomark™: Designed specifically for 3-D Tomosynthesis scans and diagnosis. They provide clear and accurate visualizations with the least potential for artifacts compared to other skin markers.



*Suremark® is proud to support Ready for Recovery Cancer Treatment Planners and their efforts to bring positive energy and motivation to the cancer recovery process. These fun and inspiring **"Moments That Matter"** labels are intended to add a splash of color and celebration to the challenges of cancer treatment. While our product line focuses on diagnostics and treatment, it is a blessing to also help celebrate the milestones towards recovery that turn cancer patients into cancer survivors. More information at www.cancerplanners.com.*



Liquimark

Finally, a line of fiducial markers for MRI that enhance accurate diagnosis, while keeping your budget in line. Liquimark markers are **NON-METALLIC**, **LATEX-FREE**, and perform consistently and accurately in a variety of applications.

LM-1 - 8mm round - for small field of view areas.

LM-2 - 15mm tube - for extremities, soft tissue, and breast MRI.

LM-3 - 30mm tube - for thoracic spines and larger areas of interest.

LM-4 - 10mm square liqui-pouch - for orthopedic imaging and pediatric cases.

LM-5 - 10mm x 24mm liqui-pouch - for extremities, soft tissue and breast MRI.

LM-6 - 10mm x 38mm liqui-pouch - for thoracic spines and larger areas of interest.

Packaged with 40 markers per box.

Price Code: SF



LM-1



LM-2



LM-3



LM-4



LM-5



LM-6





Tomomark™

Suremark® Tomomarks™ are semi-lucent markers specifically formulated for 3-D Breast Tomosynthesis scans and diagnosis. They provide clear and accurate visualizations with the least potential for artifact compared to other skin markers. These specially formulated **NON-METALLIC** markers perform consistently and accurately in a variety of applications.


- TOMO-1 - Nipple marker.
- TOMO-2 - Scar marker.
- TOMO-3 - Mole marker.
- TOMO-4 - Mass marker.
- TOMO-5 - Pain marker.

Packaged with 110 markers per box. Tomo-2 packaged with 336cm per box.


Tomo 5 packaged with 50 per box.

Price Code: SG (TOMO-1), SC (TOMO-2), and SD (TOMO-3),


SH (TOMO-4, & TOMO-5)



TOMO-1
Nipple Marker



TOMO-2
Scar Marker



TOMO-3
Mole Marker



TOMO-4
Mass Marker



TOMO-5
Pain Marker



VISION LINE Kids

It's hard enough for kids going through a difficult treatment plan, when they feel like they don't have a choice. Imaging procedures can be a scary experience for them, so having something lighthearted to think about can be the key to creating a positive experience for everyone involved. Vision Line™ Kids **NON-METALLIC** general use precision skin markers are **LATEX-FREE**, and come in six different cute and playful designs, allowing children to choose their favorite character.

Available in 2.5mm ball size.

Packaged with 110 markers per box.

Price Code: SG



VK-25A



VK-25B



VK-25C



VK-25D



VK-25E



VK-25F





VISION LINE™ Series

Visionmark™ CT

The Visionmark™ CT is the next generation in general use skin markers. These specially formulated **NON-METALLIC** markers perform consistently and accurately in a variety of applications. Easy Peel Tabs.

Available in 2.0mm and 2.5mm ball sizes.
Packaged with 110 markers per box.
Price Code: SG




V-20




V-25



V-30



V-40




V-50

Visionmark™ CT Plus

These specially formulated **NON-METALLIC** markers reduce read time on CT scans with virtually no artifact. These sizes and ball densities are perfect for all CT skin marking applications. Easy Peel Tabs.

Available in 3.0mm, 4.0mm and 5.0mm ball sizes.
Packaged with 50 markers per box.
Price Code: SD (V-30), SJ (V-40 and V-50)






Theramark™

Theramark™ markers stay securely where you put them throughout the entire treatment process. They are repositionable during the initial set up, giving you time to get the markers in just the right spot. The adhesive bonds stronger as time goes on. Theramarks adhere directly to immobilization masks without having to use tape. Only for use on immobilization masks. Not medical grade adhesive. Adhesive made with natural rubber latex.

Packaged with 110 markers per box.
Price Code: SG

TM-25





ComfortMark™ **D** GREAT FOR DIGITAL

After years of feedback from industry professionals and patients, the new ComfortMark™ labels have been developed to be the most conformable and comfortable skin markers ever. The flexible, stretchable material, easy-peel tab and adhesive-free center is not only comfortable on patient's skin, but the beautiful artwork is designed to help ease the stress of examinations.

Available in 1.5mm lead, or 2.5mm non-metallic ball sizes.

Packaged with 110 markers per box.

Price Code: SG






CT-23




CT-40

CT MARK™ and CT LINE™

The CT Mark™ and CT Line™ have been created specifically for CT imaging. This specially formulated **NON-METALLIC** material is the result of years of research and development, resulting in a virtually artifact-free opaque marker that is perfect for nearly all CT marking applications.



CT-W10



CT-W20

CT-Mark available in 2.3mm and 4.0mm ball sizes.
CT-23 Packaged with 110 markers per box.
CT-40 Packaged with 50 markers per box.
CT-W10 – 1.0mm line – packaged with 300 cm per box.
CT-W20 – 2.0mm line – packaged with 300 cm per box.
Price Code: SG (CT-23), SJ (CT-40), SK (CT-W10, CT-W20)



Suremark® CT-O

Like the rest of our CT products, the CT-O is opaque and artifact-free, making it perfect for indicating and marking details in CT Scans.

Packaged with 50 markers per box.
Price Code: SJ



Suremark® CT-CR3

Created by request, the CT-CR3 is another precision skin marking tool perfectly suited for CT applications.

Packaged with 50 markers per box.
Price Code: SJ



PortalMark LABELS

In lieu of permanent or temporary tattoos, apply **PortalMark**™ labels to delineate radiation treatment areas. Their premium medical grade adhesive is hypo-allergenic and latex-free, resulting in less irritation for the patient. PortalMark™ labels are constructed of waterproof, yet breathable material for maximum performance and comfort for up to two weeks. Our unique Suregrip™ applicator discs ensure easy and accurate placement of PortalMark™ labels without contaminating the adhesive.

The 200 Series is designed for use with darker skin surfaces, giving technologists a complete line of options.

PM-100 and 200 - Each sheet includes: 8 field lines, 4 corners, 2 set up points, and 2 isocenters.

PM-125 and 225 - Each sheet includes: 8 corners.

PM-150 and 250 - Each sheet includes: 18 isocenters, and 6 field lines.

All PortalMark™ labels are packaged with 50 sheets per box.


Price Code: SF



PM-100



PM-125



PM-150



PM-200



PM-225



PM-250



PM-300
(PM-400)*



PM-325
(PM-425)*



PM-350
(PM-450)*



PM-375
(PM-475)*



PM-300 Shown



PM-375 Shown

The newest members of the **PortalMark™** label family offer additional options for delineating radiation treatment areas. Their premium medical grade adhesive is hypo-allergenic and latex-free, resulting in less irritation for the patient. Constructed of waterproof, yet breathable material for maximum performance and comfort, PortalMarks hold strong for up to two weeks. Our unique Suregrip™ applicator discs ensure easy and accurate placement without contaminating the adhesive.

Packaged with 110 labels per box.

Price Code: SK

**PM-400 Series uses an alternate material and adhesion for more humid environments and particular skin surfaces.*



SL-30



SL-40



SL-50

Powermark™

This design maintains all of the advantages of our regular Suremark® labels, with larger lead ball sizes for visibility in radiation therapy.

Available in 3.0mm, 4.0mm, and 5.0mm ball sizes.

Packaged with 50 markers per box.

Price Code: SD (SL-30), SE (SL-40, SL-50)





SL-15RT



SL-20RT



SL-25RT



SL-15T



SL-20T



SL-25T

Suremark® Relief Tabs™

Mammograms are uncomfortable enough for patients without the painful removal of nipple markers. With this in mind, give your patients some relief with Suremark®'s Relief Tabs™. The unique adhesive-free center won't stick to sensitive areas of the skin.

Available in 1.5mm, 2.0mm and 2.5mm ball sizes.

Packaged with 110 markers per box.

Price Code: SG

Suremark® Tabs™

This design maintains all of the advantages of our regular Suremark® labels, with the added benefit of a "tabbed" feature. This enables natural and conforming placement on more prominent nipples or where more skin contact is preferred.

Available in 1.5mm, 2.0mm and 2.5mm ball sizes.

Packaged with 110 markers per box.

Price Code: SB



BL-15



BL-20

Suremark® BeautyLine™ Marks

This new upgraded design, maintains all of the technical advantages of our regular Suremark® labels including our famous “tab” feature, but is also produced with a soothing artistic flair to help take the edge off of stressful examinations for patients.

Available in 1.5mm, 2.0mm ball sizes.
Packaged with 110 markers per box.
Price Code: SA



SL-10



SL-12



SL-15



SL-18



SL-20



SL-25

Suremark®

The most common style for general use. Suremark® 1.0mm, 1.2mm, 1.5mm, 1.8mm, 2.0mm and 2.5mm labels are perfect for distinguishing between a nipple shadow and lesion, thereby reducing costly repeat examinations.

Available in 1.0mm, 1.2mm, 1.5mm, 1.8mm, 2.0mm, and 2.5mm ball sizes.
Packaged with 110 markers per box.
Price Code SA

Suremark® Clearmarks™ 

The **metal-free** semi-lucent markers are made of materials specially designed to clearly mark moles, entry points, and mass locations, yet light enough to enable clear visibility of all underlying details.


Available as SL-MO, SL-MOT, SL-SQ, SL-CR2, SL-T1, SL-A2.

Packaged with 110 markers per box.

Price Code: SC



SL-MO



Our BEST selling mole marker. SL-MO's semi-lucent design allows tissue directly under the marker to be clearly seen.



SL-MOT - Tabbed Mole Marker



SL-SQ



SL-CR2



SL-T1



SL-A2



Suremark® Mole Markers™

These markers are uniquely designed to locate raised moles and other skin nevi without overshadowing micro-calcifications. The radiolucent ring, when placed around a protuberance, prevents flattening due to compression.

Available with 2 ref. points (SL-M2) or 3 ref. points (SL-M3).

Packaged with 50 markers per box.

Price Code: SE US Patent Number 6,269,148



SL-M2

SL-M3

Suremark® BeautyLine™ Mark

With the same precision design as Clearmarks™, this **metal-free** semi-lucent marker also has beautiful artwork to help relieve the stress of exams.

Packaged with 110 markers per box.

Price Code: SC



BL-MO

CLEARLINE™ Semi-lucent Scar Markers

Clearline™ low-density, non-metallic, semi-lucent scar markers provide distinct outlining of scars and other areas of concern.

CL-8 – 0.8mm semi-lucent line – packaged with 336 cm per box.

CL-10 – 1.0mm semi-lucent line – packaged with 336 cm per box.


CL-15 – 1.5mm semi-lucent line – packaged with 320 cm per box.

CL-20 – 2.0mm semi-lucent line – packaged with 320 cm per box.

CL-25 – 2.5mm semi-lucent line – packaged with 203 cm per box.

CL-30 – 3.0mm semi-lucent line – packaged with 165 cm per box.


Price Code: SD



CL-8

CL-10


CL-15



CL-20

CL-25

CL-30




SL-W3
(.3mm wire)

SL-W4
(.4mm wire)

SL-W8
(.8mm wire)

SL-W12
(1.2mm wire)



BL-W4 (.4mm wire)

Suremark® Wire

This valuable tool is used for marking scar tissue, outlining palpable nodules, and opaque outlining of field lines and match lines.

Available in .3mm, and .4mm.
Packaged with 336cm per box.
Price Code: SC


Available .8mm, and 1.2mm wire sizes.
Packaged with 336cm per box.
Price Code: SE

Suremark® BeautyLine Wire

BeautyLine is the next generation of precision skin marking tools, providing precision accuracy for technologists and beautiful butterfly art for the patients, helping to ease their stress.

Available in .4mm.
Packaged with 336cm per box.
Price Code: SC





With the Scaneez™
 “Stick ‘n Scan” technology,
 you save time and money.
 Easy on, easy off.

Scaneez™

“Stick ‘n Scan” Cone Beam Markers. With Scaneez™ precision denture marking systems, there is no longer a need to duplicate dentures. Scaneez™ markers save on lab bills and chair time, allowing dentists to create the treatment plan and start the case during the consultation visit.

DM-1 through DM-7: Price Code SJ


Packaged with 110 markers per box

DM-PACK: Price Code SP

All 7 boxes together at a significant discount.

DM-STARTER: Price code FL

Kit includes two of each label.



Patented & LATEX-FREE



SF-10



VF-10

Non-sterile
for external use



Actual Size: 6.4mm



VF-20



CF-23



DM-8



DentalMark™ Marking System

Designed specifically for the dental industry. For a variety of uses in facial marking, denture marking, and guided surgery.

SF-10 is ideal for use with Conventional Cephalometrics.

VF-10 is ideal for Cone Beam CT.

VF-20 is ideal for Denture Cone Beam CT.

CF-23 is a larger ball diameter while eliminating artifact in Cone Beam CT Scans.

Packaged with 110 markers per box.

Price Code: SE (SF-10), SK (VF-10), SG (VF-20), SG (CF-23)

DM-8

DM-8 DentalMarks improve CBCT flange and denture base scanning.

The DM-8 standardized size (2mm) offers real-time precision data for more accurate decision making. The DM-8 is flat, so it can be used on the underside of dentures.

Packaged with 110 markers per box.

Price Code: SJ (DM-8)





SL-Arrow

*SL-Arrow being
used to mark location
of pain in patient*

Suremark® Arrow™

Easy-peel tabs with a 1mm x 8mm lead arrow on transparent, disposable label. The best tool to use when a directional marker is preferred. Especially useful to indicate foreign bodies, points of entry, ribs, and for use in radiation therapy.

Packaged with 100 markers per box.

Price Code: SE



SL-CR1

Suremark® Cross Reference™

The SL-CR1 is a 1/2" lead cross placed on a disposable tabbed label. Excellent for marking the central axis ray.

Packaged with 50 markers per box.

Price Code: SE



SL-Right

SL-Left



BH-5

SL-DT1

Suremark® Right™ and Left™

These disposable “right” and “left” markers are useful in a variety of situations where a disposable marker applied to the cassette is preferred, such as portable procedures, and isolation cases. Each label contains a 1/2" (1.2cm) lead letter mounted on a 7/8" (2.23cm) x 1 7/8" (4.76cm) backing, with an easy-peel edge for convenient removal after the label is used. This cassette marker is not latex-free.

Packaged with 100 markers per box.
Price Code: SF

Double Sided Foam Tape

Ideal for fastening film markers to cassettes.

Packaged with 110 strips per box.
SL-DT1 Price Code: FL

Convenient Box Holder

Holds five standard boxes, or two R and L boxes.
Place on a desktop or mount to a wall or imaging machine.

BH-5 Price Code: XE



suremark[®]



nltec SINCE 1994

NL-Tec Pty Ltd

*Unit 10/16 Yampi Way
Willetton, WA 6155
Australia*

Tel: +61 8 9259 5100

E: info@nl-tec.com

W: www.nl-tec.com